



IDEAL FOR THERAPY, FITNESS AND REHABILITATION

THE QUADRANGULAR BAMUSTA CUATRO PAD IS MADE FROM 60 MM THICK CLOSED CELLULAR FOAM MATERIAL. EXCELLENT DAMPING PROPERTIES, ALSO WHEN STRAINED SELECTIVELY. NON-SLIP DUE TO WAFER STRUCTURE, WATER RESISTENT AND EASY TO CLEAN.

AS SOON AS THE USER STEPS ON THE BAMUSTA CUATRO PAD ONE AUTOMATICALLY STARTS PRACTICING THE BALANCE. THE BODY WEIGHT CREATES A COMPRESSION WHICH CHALLENGES THE USER'S BALANCE AND COORDINATION.





48 x 39 x 6,0 cm



STARTING POSITION STAND ON THE CENTRE OF THE PAD ON ONE LEG.

WITH STANDING LEG SLIGHTLY BENT. BALANCE

THE BODY AND ENGAGE STOMACH MUSCLES.

EXERCISE BEND LEG AT KNEE AND LIFT AND LOWER AGAIN.

VARIANT LIFT OUTSTRETCHED LEG TO THE SIDE AND BACK

AND LOWER AGAIN.





UPPER THIGH, BUTTOCK AND DEEP MUSCLES

STARTING POSITION STAND ON THE PAD WITH FEET HIP-WIDTH APART AND

BALANCE BODY. FLEX WEIGHT ON HEELS AND STOMACH.

EXERCISE PULL BUTTOCKS BACK AND DOWNWARDS, BRINGING

THE UPPER BODY FORWARD.





UPPER THIGH, BUTTOCK AND DEEP MUSCLES

STARTING POSITION

BALANCE BODY. FLEX WEIGHT ON HEELS AND STOMACH.
PULL BUTTOCKS DOWNWARDS WHILE BRINGING THE
UPPER BODY TO THE FRONT. ALTERNATELY LIFT OUTS

STAND ON THE PAD WITH FEET HIP-WIDTH APART AND

TRETCHED ARMS OVER THE HEAD.

EXERCISE





INSIDE OF LEG AND ADDUCTORS

STARTING POSITION

STAND WITH FEET WIDER THAN HIP-WIDTH APART,

PLACING ONE FOOT ON THE PAD.

EXERCISE

DRAW HEELS TO MIDDLE, WITH TOES POINTING OUT WARDS AND WEIGHT ON THE HEEL. BEND LEGS AS FAR AS POSSIBLE WHILE KEEPING UPPER BODY STRAIGHT.



ERECTOR SPINAE AND DEEP MUSCLES

STARTING POSITION EXERCISE

STAND ON ONE LEG ON THE CENTRE OF THE PAD.

STRETCH OUT ARM AND LEG DIAGONALLY AND BALANCE
WEIGHT, BRINGING UPPER BODY FORWARD AS FAR AS
POSSIBLE.

TRENDY SPORT®



HIP FLEXOR, UPPER THIGH AND BUTTOCK MUSCLES

STARTING POSITION PLACE FRONT LEG ON THE PAD, BACK LEG ONE STEP

BACK HIP-WIDTH APART. BALANCE WEIGHT ON FRONT LEG, WITH THE FRONT KNEE ABOVE THE FOOT JOINT.

LOWER BACK KNEE TO THE FLOOR.

EXERCISE



LATERAL ABDOMINAL MUSCLES AND OUTSIDE OF LEG

STARTING POSITION STAND ON ONE LEG ON THE CENTRE OF THE PAD AND

BALANCE BODY.

EXERCISE FIRST LIFT BENT LEG AND THEN STRETCH OUT SLOWLY

AS HIGH AS POSSIBLE, WHILE BENDING THE UPPER BODY

FORWARD.

TRENDY SPORT®



LATERAL ABDOMINAL MUSCLES

STARTING POSITION EXERCISE

MITTIG EINEN FUSS AUF DEM PAD AUSBALANCIEREN. LIFT UP THE OTHER LEG TO THE SIDE AND BEND THE OUTSTRECHED ARM WITH THE ELBOW FROM ABOVE TO THE KNEE.

MAX. BELASTBARKEIT MAX. WEIGHT CARRYING CAPACITY

300 KG







NEUES MATERIAL!

NEW MATERIAL!

EXTREM ROBUST
EXTREMLY STURDY!



TRENDY SPORT®

EXERCISE INSTRUCTIONS FOR DOWNLOAD

YOU CAN DOWNLOAD ALL EXERCISE INSTRUCTIONS. JUST VISIT **WWW.TRENDY-SPORT.EU** AND GET YOUR EXERCISE INSTRUCTIONS WITH A "CLICK".

TRENDY SPORTS HOPES YOU HAVE FUN AND SUCCESS WHILE PRACTICING!

FIT WITH TRENDY SPORT!

TRENDY SPORT OFFERS EVERYTHING FOR FITNESS, THERAPY AND REHABILITATION. FIND THE BEST EXERCISE EQUIPMENT AND ACCESSORIES!

ARE YOU LOOKING FOR THE RIGHT FLOOR FOR YOUR SPORTS CLUB, YOUR YOGA ROOM OR YOUR FITNESS CENTER? TRENDY SPORT OFFERS A VARIETY OF SPORT FLOORS IN CHIC DESIGN.

TRENDY SPORT GMBH & CO. KG AM ERBERSLOHE 15-17 D-91154 ROTH-ECKERSMÜHLEN

FON +49 9171 9756-30 FAX +49 9171 9756-35

INFO@TRENDY-SPORT.COM WWW.TRENDY-SPORT.EU

